

GLOBUS ACTIVA 700

4 CHANNELS ELECTROTHERAPY | 271 PROGRAMS

PRODUCT FOR BEAUTY AND FITNESS, DESIGNED WITH SPECIFIC PROGRAMS FOR FIRMING, TONING AND REJUVENATION. THE SPECIAL HANDPIECE G-TRODE (SUPPLIED) MAKES ACTIVA 700 A TOP OF THE RANGE OF THE ELECTROSTIMULATION LINE.

60 BEAUTY | 58 FITNESS | 53 SPORT | 36 SERIAL SEQUENTIAL STIMULATION | 21 SPECIALS SPORT | 13 SKIN-FACE
12 MICROCURRENTS | 12 TENS | 3 REHAB | 3 INCONTINENCE

TECHNICAL FEATURES

Display

Visible area size 2,6"

Channels

4 independent (8 electrodes)

Frequency

0,3-150 Hz *

Pulse amplitude

40-450 μ s *

Power

0-120 mA per channel

Power supply

Rechargeable batteries

Size

100x160x35 mm


Weight

454 gr

* According to the type of current

EQUIPMENT

- 1 Carrying case
- 1 Stimulator ACTIVA 700
- 1 Handpiece G-Trode
- 4 Cables for electrode connection
- 2 Cables for microcurrents and ionophoresis
- 4 Self-adhesive square electrodes
- 4 Self-adhesive rectangular electrodes
- 1 Gel flask
- 1 Charger
- 1 Operating manual



CE
0476



TENS TENS

MENS MENS

3S 3S

STIM TIME RUN TIME

Last 10 LAST 10

Favorites FAVORITES

MULTI USER MULTI USER

PRO PROGRAMMABLE

STIM LOCK STIM LOCK

2+2 FUCTION 2+2

MY TRAINER MY TRAINER

AUTO STIM AUTO STIM

SyS SYNCRO STIM

Work Time WORK TIME

RECHARGEABLE BATTERIES

GLOBUS
ITALIAN EXCELLENCE

ACTIVA 700

4 CHANNELS ELECTROTHERAPY | 271 PROGRAMS

GLOBUS PROGRAMS LIST

Rev. En_10_2019

FITNESS

Firming

Bio-Pulse firming

Sculpting

Bio-Pulse sculpting

Toning

Mass building

Body sculpting

Definition

Jogging

Anaerobic fitness

Aerobic fitness

Cramp prevention

BEAUTY

Drainage

Bio-Pulse drainage

Lipolysis

Post-pregnancy lipolysis

Bio-Pulse relaxing massage

Toning massage

Energizing massage

Connective massage

Swollen arms

Face capillaries

Definition

Slim figure

Skin tone improvement

Post-pregnancy drainage

Post-pregnancy firming

Breast firming

Breast sculpting

SPORT

Capillarization

Warm-up

Pre-competition warm-up

Active recovery

Maximum strength

Endurance strength

Explosive strength

Reactivity

Post-competition recovery

Hypertrophy

SPECIAL SPORT

Running

Skiing

Golf

REHAB

Quadriceps atrophy (with knee prosthesis)

Recovery after ACL surgery

Shoulder subluxation prevention

TENS and PAIN

Conventional antalgic tens

Modulated antalgic tens

Endorphinic tens

Menstrual pain

Cronich pain

Cervical pain

Muscle pain

Knee pain

Scapulohumeral periarthritis

Chronic lumbago

Bursitis-tendinitis

Osteoarthritis

G-PULSE

Breast microlifting

Cleavage microlifting

Face microlifting

Skin elasticity

Cellulite

BioSkin collagen

Wrinkles

Eye wrinkles

Expression wrinkles

Stretch marks

Nail strenghtening

Bio peeling

Hematoma

MICROCURRENTS

Epicondylitis

Scapulohumeral periarthritis

Contusion

Sciatica

Lumbago

Acute pain

Articular pain

Stiff neck

Whiplash

Shoulder sprain

Knee sprain

Tendon inflammation

INCONTINENCE

Mixed incontinence

Stress incontinence

URGE incontinence

SERIAL SEQUENTIAL STIMULATION

The "3S" program list includes 36 parameter combinations.

The "3S" programs are characterized by a delay activation of the channels 3 and 4 compared with the channels 1 and 2.

The Serial Sequential Stimulation permits to stimulate the musculature in kinetic chain thanks to the differentiated activation times of the muscular groups involved.